

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: FL 31 - ELOTEX - DOLDER

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 08 May 2014 18:53:10

Date: 08 May 2014 18:37:38

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 16.09

CE Diameter Maximum (μm): 52.56

CE Diameter D[v, 0.5]: 40.15

CE Diameter D[4,3] (μm): 35.96

CE Diameter D[v, 0.9]: 51.86


CE Diameter D[3,2] (μm): 27.1

CE Diameter STDV (μm): 7.24

Particles Counted: 542

CE Diameter RSD (%): 139.69

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com