

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Microfibre 100

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morphologi

Edited On: 23 April 2014 17:38:02

Date: 23 April 2014 17:38:02

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 2.18

CE Diameter Maximum (μm): 41.24

CE Diameter D[v, 0.5]: 8.83

CE Diameter D[4,3] (μm): 14.03

CE Diameter D[v, 0.9]: 40.84


CE Diameter D[3,2] (μm): 5.411

CE Diameter STDV (μm): 1.19

Particles Counted: 33684

CE Diameter RSD (%): 82.38

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com