CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Filler LEIDI

SOP Name: MATERIE PRIME.vsop **Edited:** False

User Name: Morphologi Edited On: 09 June 2015 18:13:52

Date: 09 June 2015 18:13:52

Analysis

SOP Analysis: 3.0 SOP Fill Holes: True

SOP Segmentation: None **SOP Trash Size**: 10

SOP Optic(s) used: 50x Particles Counted: 50833


Result

CE Diameter Minimum (µm): 0.22 CE Diameter D[n, 0.1] (µm): 0.52

CE Diameter Maximum (µm): 65.65 CE Diameter D[n, 0.5] (µm): 4.52

CE Diameter Mean (µm): 5.49 CE Diameter D[n, 0.9] (µm): 10.71

CE Diameter STDV (µm): 4.55 CE Diameter RSD (%): 82.90



Comments:

Morphologi Malvern Instruments Worldwide www.malvern.com