

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: VINNAPAS LL 3030 H (WACKER)

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 26 May 2014 10:32:11

Date: 26 May 2014 10:26:25

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[v, 0.1]: 14.9

CE Diameter Maximum (μm): 92.53

CE Diameter D[v, 0.5]: 43.16

CE Diameter D[4,3] (μm): 45.04

CE Diameter D[v, 0.9]: 71.59


CE Diameter D[3,2] (μm): 28.86

CE Diameter STDV (μm): 6.14

Particles Counted: 7647

CE Diameter RSD (%): 129.04

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 140: VINNAPAS LL 3030 H (WACKER)

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com