## CE Diameter Report - Number Distribution



Measurement Details

Sample Name: OFL 190314 160W

**SOP Name:** MATERIE PRIME.vsop **Edited:** False

**User Name: Edited On:** 07 May 2014 11:38:07

Date: 07 May 2014 11:38:07

Analysis

SOP Analysis: 3.0 SOP Fill Holes: True

**SOP Segmentation**: None **SOP Trash Size**: 10

SOP Optic(s) used: 10x Particles Counted: 11490


Result

CE Diameter Minimum (µm): 1.08 CE Diameter D[n, 0.1] (µm): 2.01

CE Diameter Maximum (µm): 203.46 CE Diameter D[n, 0.5] (µm): 4.22

CE Diameter Mean (µm): 7.23 CE Diameter D[n, 0.9] (µm): 8.99

CE Diameter STDV (µm): 16.71 CE Diameter RSD (%): 231.26



Comments:

Morphologi Malvern Instruments Worldwide www.malvern.com