

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Filler DIMANNO

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morphologi

Edited On: 28 April 2015 10:18:08

Date: 28 April 2015 10:18:08

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 5.838

CE Diameter Maximum (μm): 49.93

CE Diameter D[v, 0.5]: 13.03

CE Diameter D[4,3] (μm): 17.55

CE Diameter D[v, 0.9]: 38.53


CE Diameter D[3,2] (μm): 10.93

CE Diameter STDV (μm): 3.34

Particles Counted: 19072

CE Diameter RSD (%): 81.83

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 22: Filler DIMANNO

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com