

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Aggregato da riciclo basalto

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 26 September 2014 15:00:21

Date: 04 September 2014 15:22:25

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 2.5x

Result

CE Diameter Minimum (μm): 4.29

CE Diameter D[v, 0.1]: 53.94

CE Diameter Maximum (μm): 1355.26

CE Diameter D[v, 0.5]: 492

CE Diameter D[4,3] (μm): 597.6

CE Diameter D[v, 0.9]: 1337


CE Diameter D[3,2] (μm): 126.6

CE Diameter STDV (μm): 11.21

Particles Counted: 521729

CE Diameter RSD (%): 101.65

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com