

# CE Diameter Report - Number Distribution



## Measurement Details

**Sample Name:** VINNAPAS 8031 H (WACKER)

**SOP Name:** MATERIE PRIME.v SOP

**Edited:** False

**User Name:**

**Edited On:** 28 May 2014 10:33:41

**Date:** 28 May 2014 10:33:41

## Analysis

**SOP Analysis:** 3.0

**SOP Fill Holes:** True

**SOP Segmentation:** None

**SOP Trash Size:** 10

**SOP Optic(s) used:** 20x

**Particles Counted:** 8658

## Result

**CE Diameter Minimum ( $\mu\text{m}$ ):** 0.54

**CE Diameter D[n, 0.1] ( $\mu\text{m}$ ):** 0.99

**CE Diameter Maximum ( $\mu\text{m}$ ):** 140.46

**CE Diameter D[n, 0.5] ( $\mu\text{m}$ ):** 3.37


**CE Diameter Mean ( $\mu\text{m}$ ):** 5.90

**CE Diameter D[n, 0.9] ( $\mu\text{m}$ ):** 12.49

**CE Diameter STDV ( $\mu\text{m}$ ):** 8.15

**CE Diameter RSD (%):** 138.15

CE Diameter ( $\mu\text{m}$ ) smoothed over 11 points



Record 161: VINNAPAS 8031 H (WACKER)

## Comments:

Morphologi

Malvern Instruments Worldwide

[www.malvern.com](http://www.malvern.com)