

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Formiato Neuwendis

SOP Name: 5X.vsop

Edited: False

User Name: Morphologi

Edited On: 28 August 2015 10:52:26

Date: 28 August 2015 10:52:26

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 5x

Result

CE Diameter Minimum (μm): 2.16

CE Diameter D[v, 0.1]: 62.97

CE Diameter Maximum (μm): 599.11

CE Diameter D[v, 0.5]: 183

CE Diameter D[4,3] (μm): 199.9

CE Diameter D[v, 0.9]: 339


CE Diameter D[3,2] (μm): 124.4

CE Diameter STDV (μm): 26.11

Particles Counted: 42917

CE Diameter RSD (%): 143.88

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com