

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Calce idraulica NHL 3,5

SOP Name: Calci.vsop

Edited: True

User Name: roberto antonio

Edited On: 19 June 2014 11:20:41

Date: 19 June 2014 11:09:40

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 6.419

CE Diameter Maximum (μm): 57.07

CE Diameter D[v, 0.5]: 14.76

CE Diameter D[4,3] (μm): 18.87

CE Diameter D[v, 0.9]: 37.96


CE Diameter D[3,2] (μm): 12.08

CE Diameter STDV (μm): 3.69

Particles Counted: 41945

CE Diameter RSD (%): 97.00

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 125: Calce idraulica NHL 3,5

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com