

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: KARTNOL 17.88

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 12 May 2014 11:40:39

Date: 12 May 2014 11:08:04

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 10x

Result

CE Diameter Minimum (μm): 1.08

CE Diameter D[v, 0.1]: 45.51

CE Diameter Maximum (μm): 206.31

CE Diameter D[v, 0.5]: 96.39

CE Diameter D[4,3] (μm): 108.3

CE Diameter D[v, 0.9]: 202.7


CE Diameter D[3,2] (μm): 75.07

CE Diameter STDV (μm): 14.05

Particles Counted: 9483

CE Diameter RSD (%): 170.90

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 17: KARTNOL 17.88

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com