CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Calce Idrata - Vesci

SOP Name: Calci.vsop Edited: True

User Name: roberto antonio Edited On: 18 June 2014 17:48:28

Date: 18 June 2014 17:29:22

Analysis

SOP Analysis: 3.0 SOP Fill Holes: True

SOP Segmentation: None **SOP Trash Size**: 10

SOP Optic(s) used: 50x Particles Counted: 40627


Result

CE Diameter Minimum (µm): 0.22 CE Diameter D[n, 0.1] (µm): 0.25

CE Diameter Maximum (µm): 59.05 CE Diameter D[n, 0.5] (µm): 2.42

CE Diameter Mean (µm): 3.59 CE Diameter D[n, 0.9] (µm): 8.03

CE Diameter STDV (µm): 3.95 CE Diameter RSD (%): 110.30



Comments:

Morphologi Malvern Instruments Worldwide www.malvern.com