

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: SLAG CEMENT

SOP Name: CEMENTO CALME 42.5 R TIPO II.vsop

Edited: False

User Name: Morphologi

Edited On: 17 January 2014 14:17:24

Date: 17 January 2014 14:17:24

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 6.761

CE Diameter Maximum (μm): 51.05

CE Diameter D[v, 0.5]: 18.5

CE Diameter D[4,3] (μm): 19.58

CE Diameter D[v, 0.9]: 33.17


CE Diameter D[3,2] (μm): 13.21

CE Diameter STDV (μm): 3.70

Particles Counted: 15761

CE Diameter RSD (%): 99.42

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 55: SLAG CEMENT

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com