

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Sodio Resinato

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name: Unknown

Edited On: 29 April 2014 11:50:59

Date: 29 April 2014 11:28:02

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 10635

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 1.12

CE Diameter Maximum (μm): 127.31

CE Diameter D[n, 0.5] (μm): 3.49


CE Diameter Mean (μm): 6.58

CE Diameter D[n, 0.9] (μm): 14.92

CE Diameter STDV (μm): 8.86

CE Diameter RSD (%): 134.62

CE Diameter (μm) smoothed over 11 points



Record 86: Sodio Resinato

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com