

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Filler Marmifera mineraria

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morphologi

Edited On: 27 April 2015 11:55:51

Date: 27 April 2015 11:55:51

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Particles Counted: 27309

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[n, 0.1] (μm): 0.30

CE Diameter Maximum (μm): 77.97


CE Diameter D[n, 0.5] (μm): 3.81

CE Diameter Mean (μm): 4.63

CE Diameter D[n, 0.9] (μm): 9.49

CE Diameter STDV (μm): 4.13

CE Diameter RSD (%): 89.26



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com