

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Bermacoll_ME_1000X[Chem]

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 10 June 2014 16:37:03

Date: 10 June 2014 16:28:46

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 83

Result

CE Diameter Minimum (μm): 10.36

CE Diameter D[n, 0.1] (μm): 12.58

CE Diameter Maximum (μm): 81.33

CE Diameter D[n, 0.5] (μm): 31.31


CE Diameter Mean (μm): 33.65

CE Diameter D[n, 0.9] (μm): 59.98

CE Diameter STDV (μm): 19.10

CE Diameter RSD (%): 56.78

CE Diameter (μm) smoothed over 11 points



Record 58: Bermacoll_ME_1000X[Chem]

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com