

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: E-Spheres SLG

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morfologi

Edited On: 24 April 2014 17:52:44

Date: 24 April 2014 17:52:44

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 10x

Result

CE Diameter Minimum (μm): 1.08

CE Diameter D[v, 0.1]: 87.92

CE Diameter Maximum (μm): 226.91

CE Diameter D[v, 0.5]: 152.4

CE Diameter D[4,3] (μm): 151.2

CE Diameter D[v, 0.9]: 214.3


CE Diameter D[3,2] (μm): 124.1

CE Diameter STDV (μm): 19.08

Particles Counted: 3132

CE Diameter RSD (%): 223.99

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 71: E-Spheres SLG

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com