

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Esacol 205

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Commisso

Edited On: 17 March 2014 14:28:57

Date: 17 March 2014 14:28:57

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 120595

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 1.11

CE Diameter Maximum (μm): 202.56

CE Diameter D[n, 0.5] (μm): 2.24


CE Diameter Mean (μm): 3.08

CE Diameter D[n, 0.9] (μm): 5.24

CE Diameter STDV (μm): 4.17

CE Diameter RSD (%): 135.55

CE Diameter (μm) smoothed over 11 points



Record 278: Esacol 205

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com