

# CE Diameter Report - Volume Distribution



## Measurement Details

**Sample Name:** Eloset 542

**SOP Name:** MATERIE PRIME.vsop

**Edited:** False

**User Name:** Morphologi

**Edited On:** 10 April 2014 20:15:38

**Date:** 10 April 2014 20:15:38

## Analysis

**SOP Analysis:** 3.0

**SOP Fill Holes:** True

**SOP Segmentation:** None

**SOP Trash Size:** 10

**SOP Optic(s) used:** 20x, 50x

## Result

**CE Diameter Minimum ( $\mu\text{m}$ ):** 0.22

**CE Diameter D[v, 0.1]:** 10.06

**CE Diameter Maximum ( $\mu\text{m}$ ):** 52.38

**CE Diameter D[v, 0.5]:** 32.83

**CE Diameter D[4,3] ( $\mu\text{m}$ ):** 32.65

**CE Diameter D[v, 0.9]:** 51.85


**CE Diameter D[3,2] ( $\mu\text{m}$ ):** 12.91

**CE Diameter STDV ( $\mu\text{m}$ ):** 1.07

**Particles Counted:** 14283

**CE Diameter RSD (%):** 105.92

Volume transformation: CE Diameter ( $\mu\text{m}$ ) smoothed over 11 points



Record 574: Eloset 542

## Comments:

Morphologi

Malvern Instruments Worldwide

[www.malvern.com](http://www.malvern.com)