

# CE Diameter Report - Number Distribution



## Measurement Details

**Sample Name:** Tylose\_MHS150003\_P4\_CMHSP4

**SOP Name:** MATERIE PRIME.vsop

**Edited:** True

**User Name:** Russo

**Edited On:** 10 March 2014 12:43:57

**Date:** 10 March 2014 12:14:26

## Analysis

**SOP Analysis:** 3.0

**SOP Fill Holes:** True

**SOP Segmentation:** None

**SOP Trash Size:** 10

**SOP Optic(s) used:** 50x

**Particles Counted:** 2679

## Result

**CE Diameter Minimum ( $\mu\text{m}$ ):** 0.22

**CE Diameter D[n, 0.1] ( $\mu\text{m}$ ):** 0.24

**CE Diameter Maximum ( $\mu\text{m}$ ):** 75.44

**CE Diameter D[n, 0.5] ( $\mu\text{m}$ ):** 2.15


**CE Diameter Mean ( $\mu\text{m}$ ):** 3.88

**CE Diameter D[n, 0.9] ( $\mu\text{m}$ ):** 7.93

**CE Diameter STDV ( $\mu\text{m}$ ):** 6.41

**CE Diameter RSD (%):** 165.51

CE Diameter ( $\mu\text{m}$ ) smoothed over 11 points



Record 133: Tylose\_MHS150003\_P4\_CMHSP4

## Comments:

Morphologi

Malvern Instruments Worldwide

[www.malvern.com](http://www.malvern.com)