

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Sipernat 22

SOP Name: 10X.vsop

Edited: False

User Name: Morphologi

Edited On: 28 September 2015 14:04:52

Date: 28 September 2015 14:04:52

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 10x

Result

CE Diameter Minimum (μm): 1.08

CE Diameter D[v, 0.1]: 50.76

CE Diameter Maximum (μm): 309.75

CE Diameter D[v, 0.5]: 112.5

CE Diameter D[4,3] (μm): 125.9

CE Diameter D[v, 0.9]: 215.9


CE Diameter D[3,2] (μm): 85.3

CE Diameter STDV (μm): 13.89

Particles Counted: 102694

CE Diameter RSD (%): 190.65

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 1: Sipernat 22

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com