

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: VINNAPAS 8034 H (WACKER)

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 28 May 2014 16:45:32

Date: 28 May 2014 16:29:53

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 27066

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 0.99

CE Diameter Maximum (μm): 99.61


CE Diameter D[n, 0.5] (μm): 3.85

CE Diameter Mean (μm): 5.42

CE Diameter D[n, 0.9] (μm): 10.91

CE Diameter STDV (μm): 5.59

CE Diameter RSD (%): 103.01



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com