

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: PENTACHEM - PENTARAPID CR 1001

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name: Unknown

Edited On: 13 May 2014 16:47:01

Date: 13 May 2014 16:34:11

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[v, 0.1]: 7.199

CE Diameter Maximum (μm): 98.28

CE Diameter D[v, 0.5]: 18.99

CE Diameter D[4,3] (μm): 27.84

CE Diameter D[v, 0.9]: 66.44


CE Diameter D[3,2] (μm): 14.63

CE Diameter STDV (μm): 4.09

Particles Counted: 99232

CE Diameter RSD (%): 84.01

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 37: PENTACHEM - PENTARAPID CR 1001

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com