

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: MP 2080

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name: Morphologi

Edited On: 11 April 2014 15:58:56

Date: 11 April 2014 14:01:45

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x, 50x

Particles Counted: 8872

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[n, 0.1] (μm): 0.24

CE Diameter Maximum (μm): 53.59

CE Diameter D[n, 0.5] (μm): 0.82


CE Diameter Mean (μm): 0.71

CE Diameter D[n, 0.9] (μm): 0.92

CE Diameter STDV (μm): 0.72

CE Diameter RSD (%): 101.54

CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com