## CE Diameter Report - Number Distribution



## Measurement Details

Sample Name: Calce idraulica NHL 3,5

SOP Name: Calci.vsop Edited: True

User Name: roberto antonio Edited On: 19 June 2014 11:20:41

Date: 19 June 2014 11:09:40

Analysis

SOP Analysis: 3.0 SOP Fill Holes: True

**SOP Segmentation**: None **SOP Trash Size**: 10

SOP Optic(s) used: 50x Particles Counted: 41945


Result

CE Diameter Minimum (µm): 0.22 CE Diameter D[n, 0.1] (µm): 0.46

CE Diameter Maximum (µm): 57.07 CE Diameter D[n, 0.5] (µm): 2.72

CE Diameter Mean (µm): 3.80 CE Diameter D[n, 0.9] (µm): 8.21

CE Diameter STDV (µm): 3.69 CE Diameter RSD (%): 97.00



Comments:

Morphologi Malvern Instruments Worldwide www.malvern.com