

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Mefisto L05[Chem]

SOP Name: MATERIE PRIME.vsop

User Name: roberto antonio

Date: 23 June 2014 18:11:18

Edited: False

Edited On: 23 June 2014 18:11:18

Analysis

SOP Analysis: 3.0

SOP Segmentation: None

SOP Optic(s) used: 50x

SOP Fill Holes: True

SOP Trash Size: 10

Result

CE Diameter Minimum (μm): 7.13

CE Diameter Maximum (μm): 17.28

CE Diameter D[4,3] (μm): 11.02

CE Diameter D[3,2] (μm): 10.54

Particles Counted: 715

CE Diameter D[v, 0.1]: 7.952


CE Diameter D[v, 0.5]: 10.71

CE Diameter D[v, 0.9]: 14.35

CE Diameter STDV (μm): 1.97

CE Diameter RSD (%): 20.33

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 148: Mefisto L05[Chem]

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com