

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Scagliola

SOP Name: GESSI.vsop

User Name: Antonio Roberto

Date: 16 June 2014 09:46:34

Edited: True

Edited On: 16 June 2014 09:52:30

Analysis

SOP Analysis: 3.0

SOP Segmentation: None

SOP Optic(s) used: 50x

SOP Fill Holes: True

SOP Trash Size: 10

Particles Counted: 6106

Result

CE Diameter Minimum (μm): 0.22

CE Diameter Maximum (μm): 58.35

CE Diameter Mean (μm): 4.69


CE Diameter STDV (μm): 4.59

CE Diameter D[n, 0.1] (μm): 0.67

CE Diameter D[n, 0.5] (μm): 3.52

CE Diameter D[n, 0.9] (μm): 9.49

CE Diameter RSD (%): 97.88



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com