

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: ARBOCEL PWC 50-33 LD

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 21 May 2014 18:49:37

Date: 21 May 2014 18:45:49

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 2.5x

Result

CE Diameter Minimum (μm): 4.29

CE Diameter D[v, 0.1]: 94

CE Diameter Maximum (μm): 453.77

CE Diameter D[v, 0.5]: 200

CE Diameter D[4,3] (μm): 213.8

CE Diameter D[v, 0.9]: 344.2


CE Diameter D[3,2] (μm): 148.4

CE Diameter STDV (μm): 23.83

Particles Counted: 31938

CE Diameter RSD (%): 168.13

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 114: ARBOCEL PWC 50-33 LD

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com