

# CE Diameter Report - Volume Distribution



## Measurement Details

**Sample Name:** OMYACARB 2UM --

**SOP Name:** MATERIE PRIME.vsop

**Edited:** True

**User Name:**

**Edited On:** 27 May 2014 16:46:24

**Date:** 27 May 2014 16:35:33

## Analysis

**SOP Analysis:** 3.0

**SOP Fill Holes:** True

**SOP Segmentation:** None

**SOP Trash Size:** 10

**SOP Optic(s) used:** 50x

## Result

**CE Diameter Minimum ( $\mu\text{m}$ ):** 0.22

**CE Diameter D[v, 0.1]:** 3.963

**CE Diameter Maximum ( $\mu\text{m}$ ):** 29.20

**CE Diameter D[v, 0.5]:** 8.278

**CE Diameter D[4,3] ( $\mu\text{m}$ ):** 9.06

**CE Diameter D[v, 0.9]:** 14.66


**CE Diameter D[3,2] ( $\mu\text{m}$ ):** 6.896

**CE Diameter STDV ( $\mu\text{m}$ ):** 2.39

**Particles Counted:** 22660

**CE Diameter RSD (%):** 86.56

Volume transformation: CE Diameter ( $\mu\text{m}$ ) smoothed over 11 points



Record 157: OMYACARB 2UM --

## Comments:

Morphologi

Malvern Instruments Worldwide

[www.malvern.com](http://www.malvern.com)