

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: OMYACARB 2UM --

SOP Name: MATERIE PRIME.vsop

User Name:

Date: 27 May 2014 16:35:33

Edited: True

Edited On: 27 May 2014 16:46:24

Analysis

SOP Analysis: 3.0

SOP Segmentation: None

SOP Optic(s) used: 50x

SOP Fill Holes: True

SOP Trash Size: 10

Particles Counted: 22660

Result

CE Diameter Minimum (μm): 0.22

CE Diameter Maximum (μm): 29.20

CE Diameter Mean (μm): 2.76

CE Diameter STDV (μm): 2.39


CE Diameter D[n, 0.1] (μm): 0.38

CE Diameter D[n, 0.5] (μm): 2.14

CE Diameter D[n, 0.9] (μm): 5.78

CE Diameter RSD (%): 86.56

CE Diameter (μm) smoothed over 11 points



Record 157: OMYACARB 2UM --

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com