

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: FL 31 - ELOTEX - DOLDER

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 08 May 2014 18:53:10

Date: 08 May 2014 18:37:38

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Particles Counted: 542

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[n, 0.1] (μm): 0.25

CE Diameter Maximum (μm): 52.56

CE Diameter D[n, 0.5] (μm): 3.40


CE Diameter Mean (μm): 5.19

CE Diameter D[n, 0.9] (μm): 11.07

CE Diameter STDV (μm): 7.24

CE Diameter RSD (%): 139.69

CE Diameter (μm) smoothed over 11 points



Record 3: FL 31 - ELOTEX - DOLDER

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com