

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: MP 2100

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name: Morphologi

Edited On: 14 April 2014 14:33:09

Date: 14 April 2014 13:58:10

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[v, 0.1]: 25.42

CE Diameter Maximum (μm): 178.30

CE Diameter D[v, 0.5]: 55.74

CE Diameter D[4,3] (μm): 73.75

CE Diameter D[v, 0.9]: 174.6


CE Diameter D[3,2] (μm): 46.99

CE Diameter STDV (μm): 15.07

Particles Counted: 3008

CE Diameter RSD (%): 95.70

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 596: MP 2100

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com