

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Microsfere ESPHERES SL 300

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Monphologi

Edited On: 24 April 2014 14:09:12

Date: 24 April 2014 14:09:12

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 3323

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 0.78

CE Diameter Maximum (μm): 101.60


CE Diameter D[n, 0.5] (μm): 2.71

CE Diameter Mean (μm): 4.10

CE Diameter D[n, 0.9] (μm): 8.57

CE Diameter STDV (μm): 4.65

CE Diameter RSD (%): 113.35



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com