

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: ESACOL HS 26S

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name: Antonio budriesi

Edited On: 04 June 2014 11:18:33

Date: 04 June 2014 10:52:49

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 2955

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 1.10

CE Diameter Maximum (μm): 117.74


CE Diameter D[n, 0.5] (μm): 3.82

CE Diameter Mean (μm): 6.90

CE Diameter D[n, 0.9] (μm): 12.64

CE Diameter STDV (μm): 10.88

CE Diameter RSD (%): 157.55



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com