

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Bermacoll_ME_1000X[Chem]

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 10 June 2014 16:37:03

Date: 10 June 2014 16:28:46

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Result

CE Diameter Minimum (μm): 10.36

CE Diameter D[v, 0.1]: 38.94

CE Diameter Maximum (μm): 81.33

CE Diameter D[v, 0.5]: 59.63

CE Diameter D[4,3] (μm): 58.27

CE Diameter D[v, 0.9]: 74.46


CE Diameter D[3,2] (μm): 52.65

CE Diameter STDV (μm): 19.10

Particles Counted: 83

CE Diameter RSD (%): 56.78

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com