

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Di Manno Ventilato

SOP Name: 20X.vsop

Edited: False

User Name: Morphologi

Edited On: 10 June 2015 10:19:31

Date: 10 June 2015 10:19:30

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 341358

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 1.38

CE Diameter Maximum (μm): 129.52

CE Diameter D[n, 0.5] (μm): 3.49


CE Diameter Mean (μm): 4.39

CE Diameter D[n, 0.9] (μm): 7.93

CE Diameter STDV (μm): 3.72

CE Diameter RSD (%): 84.71

CE Diameter (μm) smoothed over 11 points



Record 4: Di Manno Ventilato

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com