

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: DRY D1 C NEW GENERATION

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name: Morphologi

Edited On: 29 April 2014 08:30:19

Date: 28 April 2014 20:09:44

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x, 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 5.666

CE Diameter Maximum (μm): 54.31

CE Diameter D[v, 0.5]: 16.62

CE Diameter D[4,3] (μm): 21.54

CE Diameter D[v, 0.9]: 43.74


CE Diameter D[3,2] (μm): 11.97

CE Diameter STDV (μm): 3.31

Particles Counted: 19058

CE Diameter RSD (%): 74.32

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 82: DRY D1 C NEW GENERATION

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com