

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: TITAN 8100

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: angotti

Edited On: 17 April 2014 18:11:15

Date: 17 April 2014 18:11:15

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 4506

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 1.93

CE Diameter Maximum (μm): 133.88

CE Diameter D[n, 0.5] (μm): 6.83


CE Diameter Mean (μm): 12.70

CE Diameter D[n, 0.9] (μm): 32.71

CE Diameter STDV (μm): 14.64

CE Diameter RSD (%): 115.28

CE Diameter (μm) smoothed over 11 points



Record 12: TITAN 8100

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com