

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Eloset 542

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morphologi

Edited On: 10 April 2014 20:15:38

Date: 10 April 2014 20:15:38

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x, 50x

Particles Counted: 14283

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[n, 0.1] (μm): 0.25

CE Diameter Maximum (μm): 52.38

CE Diameter D[n, 0.5] (μm): 1.01


CE Diameter Mean (μm): 1.01

CE Diameter D[n, 0.9] (μm): 1.60

CE Diameter STDV (μm): 1.07

CE Diameter RSD (%): 105.92

CE Diameter (μm) smoothed over 11 points



Record 574: Eloset 542

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com