

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Elotex FX3300[Chem]

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Russo

Edited On: 05 March 2014 12:59:02

Date: 05 March 2014 12:59:02

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 1.04

CE Diameter D[v, 0.1]: 6.59

CE Diameter Maximum (μm): 27.24

CE Diameter D[v, 0.5]: 15.66

CE Diameter D[4,3] (μm): 16.03

CE Diameter D[v, 0.9]: 26.35


CE Diameter D[3,2] (μm): 12.41

CE Diameter STDV (μm): 4.10

Particles Counted: 405

CE Diameter RSD (%): 75.64

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 89: Elotex FX3300[Chem]

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com