

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: OMYACARB 2UM

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 27 May 2014 14:26:23

Date: 27 May 2014 14:18:00

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Particles Counted: 63591

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[n, 0.1] (μm): 0.76

CE Diameter Maximum (μm): 33.57


CE Diameter D[n, 0.5] (μm): 2.79

CE Diameter Mean (μm): 3.74

CE Diameter D[n, 0.9] (μm): 7.87

CE Diameter STDV (μm): 3.29

CE Diameter RSD (%): 87.96



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com