

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: PENTACHEM PENTARAPID AG4

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 12 May 2014 17:11:37

Date: 12 May 2014 16:58:31

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 12594

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 0.63

CE Diameter Maximum (μm): 106.90

CE Diameter D[n, 0.5] (μm): 4.45


CE Diameter Mean (μm): 6.46

CE Diameter D[n, 0.9] (μm): 14.61

CE Diameter STDV (μm): 7.31

CE Diameter RSD (%): 113.05

CE Diameter (μm) smoothed over 11 points



Record 25: PENTACHEM PENTARAPID AG4

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com