

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Calce Idrata - Vesci

SOP Name: Calci.vsop

Edited: True

User Name: roberto antonio

Edited On: 18 June 2014 17:48:28

Date: 18 June 2014 17:29:22

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 6.962

CE Diameter Maximum (μm): 59.05

CE Diameter D[v, 0.5]: 17.69

CE Diameter D[4,3] (μm): 20.41

CE Diameter D[v, 0.9]: 36.67


CE Diameter D[3,2] (μm): 13.48

CE Diameter STDV (μm): 3.95

Particles Counted: 40627

CE Diameter RSD (%): 110.30

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 120: Calce Idrata - Vesci

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com