

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Sipernat 22

SOP Name: 10X.vsop

Edited: False

User Name: Morphologi

Edited On: 28 September 2015 14:04:52

Date: 28 September 2015 14:04:52

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 10x

Particles Counted: 102694

Result

CE Diameter Minimum (μm): 1.08

CE Diameter D[n, 0.1] (μm): 1.35

CE Diameter Maximum (μm): 309.75


CE Diameter D[n, 0.5] (μm): 3.44

CE Diameter Mean (μm): 7.28

CE Diameter D[n, 0.9] (μm): 13.34

CE Diameter STDV (μm): 13.89

CE Diameter RSD (%): 190.65



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com