

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: HD 1500

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morphologi

Edited On: 10 April 2014 18:02:39

Date: 10 April 2014 18:02:39

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 10591

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 0.69

CE Diameter Maximum (μm): 135.96

CE Diameter D[n, 0.5] (μm): 1.65


CE Diameter Mean (μm): 4.66

CE Diameter D[n, 0.9] (μm): 10.78

CE Diameter STDV (μm): 9.03

CE Diameter RSD (%): 193.60

CE Diameter (μm) smoothed over 11 points



Record 573: HD 1500

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com