

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: FL 3210

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morphologi

Edited On: 10 April 2014 14:08:06

Date: 10 April 2014 14:08:06

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 2995

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 0.61

CE Diameter Maximum (μm): 156.51

CE Diameter D[n, 0.5] (μm): 6.35


CE Diameter Mean (μm): 8.45

CE Diameter D[n, 0.9] (μm): 18.31

CE Diameter STDV (μm): 10.21

CE Diameter RSD (%): 120.81

CE Diameter (μm) smoothed over 11 points



Record 568: FL 3210

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com