

# CE Diameter Report - Number Distribution



## Measurement Details

**Sample Name:** FX 7000

**SOP Name:** MATERIE PRIME.vsop

**Edited:** False

**User Name:** Morfologi

**Edited On:** 17 April 2014 09:16:40

**Date:** 17 April 2014 09:16:40

## Analysis

**SOP Analysis:** 3.0

**SOP Fill Holes:** True

**SOP Segmentation:** None

**SOP Trash Size:** 10

**SOP Optic(s) used:** 20x

**Particles Counted:** 10559

## Result

**CE Diameter Minimum ( $\mu\text{m}$ ):** 0.54

**CE Diameter D[n, 0.1] ( $\mu\text{m}$ ):** 0.98

**CE Diameter Maximum ( $\mu\text{m}$ ):** 91.38


**CE Diameter D[n, 0.5] ( $\mu\text{m}$ ):** 2.27

**CE Diameter Mean ( $\mu\text{m}$ ):** 5.47

**CE Diameter D[n, 0.9] ( $\mu\text{m}$ ):** 14.21

**CE Diameter STDV ( $\mu\text{m}$ ):** 7.91

**CE Diameter RSD (%):** 144.49



**Comments:**

Morphologi

Malvern Instruments Worldwide

[www.malvern.com](http://www.malvern.com)