

# CE Diameter Report - Number Distribution



## Measurement Details

**Sample Name:** Formiato Shandong

**SOP Name:** 10X.vsop

**Edited:** False

**User Name:** Morphologi

**Edited On:** 27 November 2015 13:10:51

**Date:** 27 November 2015 13:10:51

## Analysis

**SOP Analysis:** 3.0

**SOP Fill Holes:** True

**SOP Segmentation:** None

**SOP Trash Size:** 10

**SOP Optic(s) used:** 10x

**Particles Counted:** 22677

## Result

**CE Diameter Minimum ( $\mu\text{m}$ ):** 1.08

**CE Diameter D[n, 0.1] ( $\mu\text{m}$ ):** 2.06

**CE Diameter Maximum ( $\mu\text{m}$ ):** 321.95


**CE Diameter D[n, 0.5] ( $\mu\text{m}$ ):** 9.05

**CE Diameter Mean ( $\mu\text{m}$ ):** 22.24

**CE Diameter D[n, 0.9] ( $\mu\text{m}$ ):** 61.96

**CE Diameter STDV ( $\mu\text{m}$ ):** 28.73

**CE Diameter RSD (%):** 129.17



## Comments:

Morphologi

Malvern Instruments Worldwide

[www.malvern.com](http://www.malvern.com)