

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Bactiblock 101

SOP Name: MATERIE PRIME.vsop

User Name: morfologi

Date: 22 April 2014 18:34:02

Edited: True

Edited On: 22 April 2014 18:41:25

Analysis

SOP Analysis: 3.0

SOP Segmentation: None

SOP Optic(s) used: 10x, 50x

SOP Fill Holes: True

SOP Trash Size: 10

Particles Counted: 37661

Result

CE Diameter Minimum (μm): 0.22

CE Diameter Maximum (μm): 90.60

CE Diameter Mean (μm): 2.25

CE Diameter STDV (μm): 2.11


CE Diameter D[n, 0.1] (μm): 0.62

CE Diameter D[n, 0.5] (μm): 1.72

CE Diameter D[n, 0.9] (μm): 4.12

CE Diameter RSD (%): 93.66

CE Diameter (μm) smoothed over 11 points



Record 31: Bactiblock 101

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com