

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: VINNAPAS 4410 L (Wacker)

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 23 May 2014 10:57:16

Date: 23 May 2014 10:50:19

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 10x

Result

CE Diameter Minimum (μm): 1.08

CE Diameter D[v, 0.1]: 29.08

CE Diameter Maximum (μm): 187.66

CE Diameter D[v, 0.5]: 76.79

CE Diameter D[4,3] (μm): 88.37

CE Diameter D[v, 0.9]: 164.8


CE Diameter D[3,2] (μm): 54.72

CE Diameter STDV (μm): 12.79

Particles Counted: 17637

CE Diameter RSD (%): 114.00

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com